# Simply a good feeling. Systematic hygiene from Dürr Dental.

4 colours, 4 areas - simply logical, simply safe





# It is best when you can rely on a system

Back in 1965 Dürr Dental already developed the first disinfectant for suction systems in the form of Orotol<sup>®</sup>. This marked the beginning of the continuous development and refinement of our expertise in the hygiene area. Today, dentists throughout the world trust Dürr System Hygiene to ensure comprehensive protection for their surgeries.

# Dürr System Hygiene

# Consulting

- Training courses, staff training sessions and on-site consulting
- Expert knowledge from qualified trainers
- Solution-oriented practice tips tailored to individual dental surgeries
- Support for development of a quality management system



## Service

- Professional and solution-oriented product advice
- Up-to-date online information





## Quality

- High production standards Made in Germany
- Defined and controlled work processes ensure safety
- Decades of experience in the manufacture of disinfectants and cleaning agents



In addition to Dürr System Hygiene products, we also offer a comprehensive service package to support your practice team in their daily work. We stand by your team with our on-site consulting service. Our service department will help you with questions about our products. And outstanding quality as the product of decades of experience goes hand-in-hand with frequent innovations to meet the needs of dental practices.

## User orientation

- Guidance during application thanks to the Dürr System Hygiene colour system
- Distinctive and easy-to-store 2.5 | bottle from Dürr System Hygiene
- Ergonomic bottle shape makes for easy handling
- Application aids for perfect dosing

# Hygiene technology

- Over 75 years of experience in medical technology
- Unbeatable: perfect technical solutions combined with the cleaning and disinfecting solutions from Dürr Dental
- Gentle and effective products shine with easy handling and speed of use





# It is best when a single system offers a complete solution

When it comes to infection protection, disinfection and cleaning, nothing should be left to chance. This is why we have carried out countless tests to demonstrate that all our preparations meet the key demands of the day-today surgery environment.

# Safe

- Comprehensive spectrum of activity:
- bactericidal, fungicidal, tuberculocidal, virucidal
- Certified, reliable effect
- VAH-listed and EN-compliant, without exception
- High disinfecting and cleaning power



## Gentle

- Superior combination of active agents
- Particularly gentle and kind to skin
- Suitable for use on most surfaces and recommended by leading equipment manufacturers
- Clinically tested
- Bio-degradable

## Convenient

- Easy and efficient to use
- Convenient, standard container sizes
- Cleaning according to plan
- Dosing in caps
- Highly cost-effective to use

- Fast
- Short reaction times
- Quick, residue-free drying
- Practical application Simple dosing



# Dürr Dental Hygiene technology – for maximum efficiency in the hygiene cycle

For more than 75 years, Dürr Dental has been synonymous with trend-setting developments in medical technology – made in Germany. With this tradition and the goal of systematically perfecting our range of hygiene technology products, today we can offer a variety of innovative

Use

Manual cleaning and disinfection Machine reprocessing in the washer-disinfector

Function test/check and instrument care

Packaging of parts for sterilisation with Hygopac View and Hygofol



technical solutions so that day-to-day hygiene and preparatory measures in dental surgeries can be performed most efficiently and with the best results. In Dürr System Hygiene, the individual components work together to create an unbeatable team.



# Instruments



# Disinfection and cleaning Instruments and drills

■ ID 212 Instrument disinfection Concentrate for cleaning and disinfecting instruments with the best material compatibility	D 212 Market	<ul> <li>Bactericidal incl. TB, yeasticidal, limited virucidal effect &amp; effect on non-enveloped adenoviruses and noroviruses</li> <li>Suitable for general, surgical and rotary instruments</li> <li>Very short reaction time in an immersion bath</li> <li>Excellent material compatibility with instruments containing rubber and silicon</li> <li>Active ingredients: Quaternary ammonium compounds</li> <li>Available as: 2.5 1</li> <li>Our recommendation: Hygosonic, Hygobox</li> </ul>
<b>ID 212 forte Instrument disinfection</b> Concentrate for instrument cleaning and disinfection with a short reaction time in an ultrasonic bath	D222 w Morks days D222 w Morks days days D222 w Morks days	<ul> <li>Bactericidal incl. TB, yeasticidal, limited virucidal effect &amp; effect on non-enveloped adenoviruses and noroviruses</li> <li>Suitable for general, surgical and rotary instruments</li> <li>High corrosion protection and long working time of 7 days</li> <li>Extremely fast disinfection and cleaning in an ultrasonic bath</li> <li>Active ingredients: Alkylamines, quaternary ammonium compounds</li> <li>Available as: 2.5 1</li> <li>Our recommendation: Hygosonic, Hygobox</li> </ul>
■ ID 213 Instrument disinfection Fully virucidal concentrate for cleaning and disinfecting general, surgical and rotating instruments	P 23 P 23	<ul> <li>Bactericidal (incl. TB), fungicidal, virucidal (certified to EN 14476)</li> <li>Exceptional cleaning efficiency - dissolves even the most stubborn contamination</li> <li>Very economical thanks to a long working time of 14 days</li> <li>Corrosion inhibitors make for excellent material compatibility</li> <li>Active ingredients: Alkylamines, quaternary ammonium compounds</li> <li>Available as: 2.5 l</li> <li>Our recommendation: Hygosonic, Hygobox</li> </ul>
<b>ID 215 Enzymatic instrument cleaner</b> Concentrate for manual, enzymatic instrument cleaning prior to disinfection		<ul> <li>Bactericidal, yeasticidal, limited viricidal effect</li> <li>The microbicidal effect protects personnel</li> <li>Exceptional cleaning efficiency - dissolves even the most stubborn contamination</li> <li>Excellent material compatibility</li> <li>Active ingredients: Quaternary ammonium compounds, guanidine compounds</li> <li>Available as: 2.5  </li> <li>Our recommendation: Hygosonic, Hygobox, Bur bath</li> </ul>
■ ID 220 Bur disinfection Fully virucidal, ready-to-use solution for cleaning and disinfecting rotating instruments	E 220 W and a construction of the constructio	<ul> <li>Bactericidal (inc. TB), fungicidal, virucidal (certified to DVV/RKI and EN 14476)</li> <li>High cleaning power, especially designed for the microstructures in drills and burs</li> <li>High corrosion protection and long working time of 7 days</li> <li>Very short reaction times</li> <li>Active ingredients: Alcohols in alkaline medium</li> <li>Available as: 2.5 l</li> <li>Our recommendation: Hygosonic, Bur bath</li> </ul>



# Devices and accessories

with DIN EN ISO 11607-2.

Hygobox Sturdy and robust transport and disinfection box with a 3 litre capacity	<ul> <li>Ergonomic handles for safe transport</li> <li>Tilt-resistant drip-off position for the integrated cartridge</li> <li>Specially shaped slits prevent the instruments from becoming entangled</li> <li>Dimensions (HxWxD): 13 x 32.5 x 21 cm</li> <li>Colours: Anthracite-coloured lid with blue or yellow cartridge, lid available in 5 different colours and with a white cartridge</li> <li>Our recommendation: ID 212, ID 212 forte, ID 213, ID 215, MD 520, MD 535</li> </ul>
ID 220 bur box	<ul> <li>with lid and cartridge</li> </ul>
	made of robust plastic
Autoclavable disinfection box for rotating	<ul> <li>Autoclavable up to 121°C</li> </ul>
instruments	Dimensions: Height 8 cm, Ø: 8.8 cm
	Capacity: 200 ml
Hygosonic         Heatable ultrasonic device for thorough cleaning of instruments, aids, prostheses, etc. The instruments can also be reliably and safely disinfected using the ID products.	<ul> <li>Effective and material-protecting disinfection and cleaning</li> <li>Sound field optimisation ensures even cleaning</li> <li>Enhanced performance with 'boost function'</li> <li>De-gasification of solutions with 'de-gas function'</li> <li>Dimensions (HxWxD): 21.4 x 30 x 17.9 cm</li> <li>Trough volume: 2.75 I</li> <li>Accessories: Lid, glass set, insert / immersion basket</li> <li>Our recommendation: ID 212, ID 212 forte, ID 213, ID 215, ID 220, MD 530, MD 533</li> </ul>
	<ul> <li>Optimal pressing pressure for consistent sealing seam quality</li> </ul>
Нудорас	<ul> <li>Optimal pressing pressure for consistent sealing seam quality</li> <li>Assembly-line principle for fast and economical operation</li> </ul>
Rotary sealing device for transparent	<ul> <li>Continuous temperature regulation for optimal seal seam</li> </ul>
sterile packaging	<b>Dimensions (HxWxD):</b> 16.0 x 43.5 x 13.5 cm
	Our recommendation: Hygofol bags, Hygofol rolls
Hygopac View	<ul> <li>Exceptionally compact design with minimal space requirements</li> <li>Intuitive user interface via touch display</li> </ul>
	<ul> <li>Process can be documented by saving to an SD card and via the network</li> </ul>
Validatable rotary sealing device	connection
for standard-compliant sealing of	
sterilised packaging in accordance	<b>Dimensions (HxWxD):</b> 17.0 x 37.0 x 14.0 cm
	Our recommendation: Hyaofol baas, Hyaofol rolls, instrument table, Hyaofol station

Our recommendation: Hygofol bags, Hygofol rolls, instrument table, Hygofol station



# Devices and accessories

## Hygoclave 90

Class B steam steriliser in accordance with EN 13060 with innovative DuraSteam technology. The key features of the Hygoclave 90 are the reliability and impressively fast cycle times combined with high energy efficiency. This is achieved by the perfect combination of and trouble-free interplay between a steam condenser, three-piston vacuum pump, heat pipe, welded valve seats and powerful heating output of 2400 W.



#### Hygotest

The Helix test body with 250 test strips for a regular, quick and simple Helix test.



- Added to this are exceptional user-friendliness and ease of maintenance, which – together with the fully stainless steel design designed to maximise value retention of the instrument – help set the standards for a new era in performance.
- One of the fastest sterilisation cycles (universal program 134°C/5.5 minutes from 12 minutes)
- Low water and energy consumption
- High loading capacity (up to 10 kg when using sterile containers)
- High sterilisation safety, thanks to a strong and quickly established vacuum (high condensation rate)
- Approval incl. all relevant parameters directly at the device by signature
- Logs can be directly saved as PDFs without a PC connection
- Connection to Tyscor Pulse and third party suppliers
- Easy and fast maintenance
- The carrier can be loaded with up to 6 trays or rotated by 90° -2 standard tray cassettes and 1 sterile container (height 65 mm)

Dimensions (HxWxD): 47.5 x 54 x 52 cm

- In accordance with EN 867-5 and ISO 11140-1
- PCD (Process Challenge Device) and 250 indicators
- Can be used as a batch control system and for a Bowie-Dick test
- Starter set with 50 indicators included with the initial autoclave delivery

### Hygofol

Hygofol sterilised packaging offers maximum reliability with its multi-layer laminated foil and highly tear-resistant sterilisation paper



- High puncture resistance and excellent mechanical strength
- High sealing seam tensile strength
- Hygofol as practical, prefabricated bags, available in 3 sizes
- Hygofol in rolls, available in 5 sizes for individual packaging

#### Dimensions:

Hygofol Bag (WxL): 7.5 cm, 10 cm or 15 cm x 30 cm each Hygofol Rolls (WxL): 5 cm, 7.5 cm, 10 cm, 15 cm or 25 cm x 100 m each

\* Tested with sterile containers made of aluminium, each with two permanent filters (one in the lid and one in the base tray). Weight includes sterile containers and accompanying cartridges.







### Hygoprint Plus

Automatic label printer for sterile product packaging permits direct printing of program-specific labels.



### Flexible positioning thanks to Ethernet connection

- Label size in accordance with guidelines can be applied outside the sealing seam
- Bar code for scanning in digital data management
- Label specifications: device (REF and serial number), program, LOT date, user (approval by), expiration date

## **Dimensions (HxWxD):** 17.5 x 12.5 x 24 cm

Our recommendation: Hygoprint Plus labels and toner rolls

- Quick and simple documentation of the sterilisation process
- Complete labelling with all relevant information
- Very lightweight and handy

**Our recommendation:** Single or double-sided adhesive labels, documentation block with 100 sheets, replacement ink roller

System for labelling sterile product packaging.



## Hygodem 90

Hygoprint

Hygodem 90 delivers high quality in a compact unit. The ion exchanger unit filters tap water in three cartridges and can be directly connected to the Hygoclave 90 even in the tiniest spaces.



- Ion exchange unit
- Convenient, thanks to 3 ion exchanger cartridges
- Simple maintenance, thanks to relief valve
- Simple installation, thanks to quick connector
- Low space requirements
- Affordable purchase price

Our recommendation: Cartridge set for Hygodem 90

#### Hygopure 90

When running at maximum performance, the Hygoclave 90 needs a powerful and reliable partner to cope with the pressures of a busy day-to-day practice environment. With its high capacity, the Hygopure 90 allows large quantities (several litres) of DI water to be taken out at once.



- Reverse osmosis unit incl. removal tap and storage tank
- Multiple filter stages + high-performance membrane with a service valve for checking effectiveness
- Quick connector makes for simple installation
- High capacity
- Multiple litres of DI water can be removed at once
- The eco module permits water savings of up to 70 %

**Our recommendation:** Cartridge set for Hygopure 90, DF50 membrane for Hygopure 90, stainless steel extraction tap for Hygopure 90, HMEC1 conductance meter, LD1 Aqua-Stop valve

# Surfaces



# Disinfection and cleaning Surface disinfection

## Bactericidal incl. TB, yeasticidal, virucidal (certified to DVV/RKI and EN 14476) FD 300 Surface disinfection Especially fast reaction time Good material compatibility Fully virucidal and alcohol-free concentrate for disinfecting surfaces of medical devices Active ingredients: Alkylamines, quaternary ammonium compounds Available as: 2.5 Our recommendation: FD 370 cleaner, FD multi wipes, FD multi wipes compact FD 312 Surface disinfection\* Bactericidal, yeasticidal, limited virucidal PLUS (enveloped viruses as well as adenoviruses and noroviruses) • Economical application: up to 250 | of ready-to-use solution at 1% Alcohol-free concentrate for disinfection Certified material compatibility of surfaces and objects Pleasant scent Active ingredients: Quaternary ammonium compounds Available as: 2.5 Our recommendation: FD 370 cleaner, FD multi wipes, FD multi wipes compact FD 370 cleaner • Effective cleaning agent for surgery equipment, devices and floors Cleaning agent for medical practices Disinfectant effect not neutralised No stickiness or discolouring of cleaned surfaces in combination with disinfectants • High yield: up to 250 l ready-to-use solution at 1% Concentrate for cleaning surfaces and devices Ingredients: Non-ionic surfactants Available as: 2.5 | Our recommendation: FD 300, FD 312 Quick-acting disinfection FD 322 Quick-acting disinfection Bactericidal incl. TB, yeasticidal, fungicidal, limited virucidal PLUS (enveloped viruses as well as adenoviruses and noroviruses) Very short reaction time Fast-acting, ready-to-use solution Dries without residue for disinfecting and cleaning surfaces

Certified material compatibility

Active ingredients: Alcohols

FD 322 available as: 750 ml for Hygowipe Plus, 1 l, 2.5 l, 10 l

Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes,

FD multi wipes compact

FD 322 perfume free available as: 2.5 l, 10 l

Our recommendation: 600 ml spray bottle, FD multi wipes, FD multi wipes compact

of medical devices and fixtures



#### FD 333 Quick-acting disinfection

Fast acting, ready-to-use solution for quick-acting disinfection of surfaces of medical devices and fixtures



### FD 333 forte Quick-acting disinfection

Fully virucidal ready-to-use solution for the quick-acting disinfection of the surfaces of medical devices and medical inventory

#### FD 366 sensitive Disinfection of sensitive surfaces

Alcohol-reduced, ready-to-use solution for disinfecting and cleaning especially sensitive surfaces of medical devices



# Care

### FD 360 Cleaning and care of vinyl upholstery

Cleaning fluid for deep cleaning and care of vinyl upholstery



- Bactericidal incl. TB, yeasticidal, fungicidal, limited virucidal PLUS (enveloped viruses as well as adenoviruses and noroviruses)
- Excellent material compatibility
- Fresh lemon scent
- Quick drying

Active ingredients: Alcohols, quaternary ammonium compounds Available as: 750 ml for Hygowipe Plus, 1 l, 2.5 l, 10 l Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact

- Bactericidal, mycobactericidal incl. TB, yeasticidal, fungicidal, virucidal (certified to DVV/RKI and EN 14476)
- Good material compatibility
- Fresh lemon scent
- Quick drying

Active ingredients: Alcohols, quaternary ammonium compounds Available as: 1 l, 2.5 l, 10 l Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact

- Bactericidal incl. TB, yeasticidal, limited virucidal effect & effect on non-enveloped noroviruses
- Quick-acting disinfection effect
- Ideal for treatment units in combination with FD 360 Cleaning and care of vinyl upholstery
- Certified material compatibility recommended by leading equipment manufacturers

Active ingredients: Alcohols, quaternary ammonium compounds Available as: 750 ml for Hygowipe Plus, 1 l, 2.5 l, 10 l Our recommendation: Spray head, 600 ml spray bottle, FD multi wipes, FD multi wipes compact

- Cleaning agent, gentle on materials, combined with a care product for long-term value-retention
- Intensive cleaning effect removes stains of all kinds, rims and discolouration
- Certified material compatibility recommended by leading equipment manufacturers
- Ideal for treatment units in combination with FD 366 sensitive

**Ingredients:** Special surfactants, silicon compounds and care formula based on avocado oil

Available as: Set with 500 ml spray bottle, 1 FD polishing wipe and 9 special sponges Our recommendation: FD 366 sensitive, FD 366 sensitive wipes, FD polishing wipe

# Surface disinfection wipes

# Alcohol-free disinfection wipes



FD 300 top wipes Surface disinfection

Fully virucidal and alcohol-free disinfection and cleaning wipes for surfaces of medical devices



- Bactericidal (incl. TB), fungicidal, virucidal (certified to EN 14476)
- Especially fast reaction time
- Easy to remove from the clip closure package
- Easy-to-grip structure for quick handling

Active ingredients: Alkylamines, quaternary ammonium compounds Available as: Package with 50 wipes Wipe size: 17.5 x 25 cm

- Bactericidal, yeasticidal, limited virucidal effect
- High yielding surface coverage of up to 8 m<sup>2</sup>
- Wipes with tear-resistant honeycomb structure
- With fresh scent

Active ingredients: Quaternary ammonium compounds Available as: Package with 15 wipes Wipe size: 22 x 42 cm

Alcohol-free disinfection and

FD 312 wet wipes\* Surface disinfection

cleaning wipes for large surfaces



# Alcohol-based disinfection wipes

 FD 322 premium wipes quick-acting disinfection,
 FD 322 top wipes quick-acting disinfection

Alcohol-based, quick-acting disinfection and cleaning wipes for smaller surfaces of medical devices



FD 333 wipes quick-acting disinfection, FD 333 top wipes quick-acting disinfection

Alcohol-based, ready-to-use disinfection and cleaning wipes for smaller surfaces of medical devices



 FD 333 forte premium wipes quick-acting disinfection,
 FD 333 forte wipes quick-acting disinfection

Fully virucidal ready-to-use cleaning and disinfection wipes for the surfaces of medical devices and medical inventory



- Bactericidal incl. TB, yeasticidal, fungicidal, limited virucidal PLUS (enveloped viruses as well as adenoviruses and noroviruses)
- Very short reaction time
- Dries without residue
- Certified material compatibility

# Active ingredients: Alcohols

FD 322 premium wipes	FD 322
Available as:	Availab
Dispenser box with 80 wipes	Package
Wipe size: 14 x 19 cm	Wipe si

FD 322 top wipes Available as: Package with 50 wipes Wipe size: 17.5 x 25 cm

- Bactericidal incl. TB, yeasticidal, fungicidal, limited virucidal PLUS (enveloped viruses as well as adenoviruses and noroviruses)
- Quick-acting disinfection effect
- Quick drying times
- Perfume-free

Active ingredients: Alcohols, quaternary ammonium compounds FD 333 wipes Available as: Dispenser box with 110 wipes

FD 333 top wipes Available as: Package with 50 wipes Wipe size: 17.5 x 25 cm

NEW

- Bactericidal, mycobactericidal incl. TB, fungicidal, virucidal (certified to DVV/RKI and EN 14476)
- Quick-acting disinfection effect
- Extra-large wipes

Wipe size: 14 x 22 cm

Fresh lemon scent

Active ingredients: Alcohols, quaternary ammonium compounds

FD 333 forte premium wipes Available as: Dispenser box with 80 wipes Wipe size: 14 x 19 cm FD 333 forte wipes Available as: Dispenser box with 100 wipes Wipe size: 20 x 30 cm

#### FD 350 Disinfection wipes

Alcohol-based, quick-acting disinfection and cleaning wipes for smaller surfaces of medical devices



- Bactericidal incl. TB, yeasticidal, fungicidal, limited virucidal PLUS (enveloped viruses as well as adenoviruses and noroviruses)
- Quick-acting disinfection effect
- Dries without residue
- 3 different fragrances: Classic, Flower, Lemon

Active ingredients: Alcohols

Available as: Dispenser box with 110 wipes Wipe size: 14 x 22 cm



 FD 366 sensitive disinfection of sensitive surfaces,
 FD 366 sensitive top wipes disinfection of sensitive surfaces

Ready-to-use disinfection and cleaning wipes for sensitive surfaces of medical devices

# Dry wipes

### FD multi wipes

Universal dry wipes, can be soaked in disinfectant and cleaning agents



FD multi wipes Winter Winter

- Bactericidal, tuberculocidal, yeasticidal, limited virucidal effect & effect on non-enveloped noroviruses
- Quick-acting disinfection effect
- Extra-large wipes
- Perfume-free

Active ingredients:

Alcohols, quaternary ammonium compounds

FD 366 sensitive wipes Available as: Dispenser box with 100 wipes

Wipe size: 20 x 30 cm

FD 366 sensitive top wipes Available as: Package with 50 wipes Wipe size: 17.5 x 25 cm

- All-purpose wipes can be used together with surface disinfectants
- Lint-free and tear-proof
- Excellent wetting; does not scratch surfaces
- Very economical

Available as: Dispenser box with 180 wipes Wipe size: 20 x 30 cm Our recommendation: FD 300, FD 312, FD 322, FD 333, FD 333 forte, FD 366 sensitive, wall bracket for FD multi wipes box

FD multi wipes compact

Universal dry wipes, can be soaked in disinfectant and cleaning agents



- All-purpose wipes can be used together with surface disinfectants
- Lint-free and tear-proof
- Excellent wetting; does not scratch surfaces
- Compact box, extra-large wipes

Available as: Dispenser box with 60 wipes Wipe size: 29 x 30 cm Our recommendation: FD 300, FD 312, FD 322, FD 333, FD 333 forte, FD 366 sensitive

## FD polishing wipe

Microfibre cloth with special fibres for gentle cleaning and polishing of sensitive surfaces



- Particularly thorough thanks to the fine fibre material
- Fast, streak-free cleaning
- Antistatic properties remove small particles of dirt and dust
- Suitable for washing at temperatures of up to 60°C

Available as: Individual wipe Wipe size: 32 x 32 cm

# Skin and hands



# Disinfection and cleaning Hand disinfection

## HD 410\* Hand disinfection

Alcoholic rub-in preparation for hygienic and surgical hand disinfection



# Hand cleaning

HD 425 Cleansing lotion

Mild washing lotion to wash skin

and hands

## Bactericidal incl. TB, yeasticidal, limited virucidal effect & effect on non-enveloped adenoviruses, noroviruses and rotaviruses

- Extremely fast: hygienic hand disinfection in 15 seconds / surgical hand disinfection in 90 seconds
- Pleasant to the skin, thanks to rehydrating effect
- Perfume-free

### Active ingredients: Alcohols

Available as: 500 ml, 500 ml for Hygocare Plus, 1 l, 2.5 l, 10 l Our recommendation: Universal dispenser S 400, L 400, L 1000, touchless T 400, T 1000, RX 5M, RX 5T, dosing pumps

- Dermatologically tested very kind to skin
- Contains skin-care and lipid replenishing substances
- Ideal combination product for alcohol-based hand disinfection HD 410

**Ingredients:** Active cleaning surfactants and skin-care components, such as allantoin **Available as:** 500 ml

Our recommendation: HD 410, HD 440

HD 435 Cleansing lotion

Mild cleansing lotion with a skin-typical pH value to retain the natural skin flora



# Care

### HD 440 Care lotion

Lotion for stressed and sensitive skin



- Intensive cleaning thanks to selected active washing components
- Contains skin care and lipid replenishing substances as well as skin moisturising factors
- Ideal combination product for alcohol-based hand disinfection HD 410

Ingredients: Active cleaning surfactants, allantoin Available as: 500 ml, 1 ml, 2.5 l Our recommendation: Universal dispenser S 400, L 400, L 1000, Touchless T 400, T 1000, RX 5M, RX 5T, dosing pumps, HD 410, HD 440

- Protects and cares for skin, prevents inflammation
- Prevents drying out of the skin
- Tried and tested ingredients keep the hands smooth
- Non-greasy quickly absorbed in the skin

**Ingredients:** Nourishing oils, witch hazel extract, allantoin **Available as:** 500 ml

**Our recommendation:** Universal dispenser S 400, L 400, Touchless T 400, RX 5M, RX 5T, dosing pump, HD 410, HD 425, HD 435





# Dispenser system

## Touchless universal preparation dispenser

Battery-operated sensor dispenser for contactless, hygienic application of hand preparations

SENSOR



### Universal dispensers

Universal wall-mounted dispenser for application of hand products



Perfect hand hygiene thanks to touchless operation

- Battery operation with an acoustic warning signal that lets you know when the battery needs to be replaced
- High standard of hygiene, thanks to an anodised aluminium housing
- Autoclavable stainless steel pump or disposable pump

Dimensions (WxHxD): T 1000/T 400: 9.2 x 33.0 x 21.0 cm Available as: T 400 for 500 ml bottle, T 1000 for 1 l bottle Our recommendation: HD 410, HD 435 in a 500 ml or 1 l bottle, HD 440 in a 500 ml bottle, dish holder, table stand, dispenser column, disposable pumps

- Extremely user-friendly
- Available with short (S) or long (L) arm lever
- Durable aluminium and stainless steel version
- Autoclavable pump or disposable pump

**Dimensions (WxHxD):** S 400/L 400: 8.2 x 24.5 x 8.2 cm, arm lever length: S: 16 cm, L: 21.5 cm/L 1000: 9.2 x 29 x 9 cm, arm lever length: 22.5 cm **Available as:** S 400 and L 400 for 500 ml bottle, L 1000 for 1 l bottle **Our recommendation:** HD 410, HD 435 in a 500 ml or 1 l bottle, HD 440 in a 500 ml bottle, disposable pumps

Universal dispenser RX

Manual or battery-powered wall-mounted universal dispenser for application of hand products



- Manual (RX 5M) or touchless (RX 5T)
- Housing and collector tray made of robust white plastic
- Easy pump replacement simply remove from the front

Dimensions (WxHxD): 10.2 x 27.0 x 14.2 cm Available as: RX 5M and RX 5T for 500 ml bottle Our recommendation: HD 410, HD 435, HD 440 in a 500 ml bottle

\*Use disinfectants with caution. Always read the label and product information before use.

# Special areas

# Disinfection and cleaning Suction system and spittoon

Orotol<sup>®</sup> plus disinfection of suction systems

Foam-free liquid concentrate for simultaneous cleaning, disinfection, deodorisation and care of all suction systems and amalgam separators

#### Orotol<sup>®</sup> ultra disinfection of suction systems

Powder concentrate for simultaneous cleaning, disinfection, deodorisation and care of all suction systems and amalgam separators

## MD 555 cleaner Special detergent for suction systems

Foam-free cleaning concentrate for all types of suction system and amalgam separator

### MD 550 Spittoon bowl cleaner

Foam-free, ready-to-use solution for intensive, residue-free cleaning of spittoons

## OroCup Care system

User-friendly, closed care system for disinfection and cleaning of suction systems and amalgam separators

# Problem-free removal of residues, such as lime scale, blood or plaque colouring

#### Active ingredients: Alcohols

Available as: 750 ml

Our recommendation: Orotol® plus, Orotol® ultra, MD 555 cleaner

- Creates an air-water mix for optimal wetting and effective cleaning of the suction system
- Efficient work one system for application and suctioning of the ready-to-use solutions
- Connecting adapter for all suction hoses

#### Available as: 2 litre care system

Our recommendation: Orotol® plus, Orotol® ultra, MD 555 cleaner

Bactericidal, yeasticidal, limited virucidal effect & effect on non-enveloped adenoviruses and noroviruses

- Dissolves biofilm and prevents formation of sludge through blood and proteins
- The solution for contaminated suction systems provides reliable protection against infection
- Certified material compatibility recommended by leading equipment manufacturers

Active ingredients: Quaternary ammonium compounds Available as: 1 |, 2.5 | Our recommendation: OroCup, MD 555 cleaner, MD 550

- Bactericidal incl. TB, sporicidal (effective against Clostridium difficile), yeasticidal, fungicidal, virucidal (certified to DVV/RKI and EN 14476)
- Dissolves biofilm and prevents formation of sludge through blood and proteins
- Low application concentration (1%)
- Certified material compatibility recommended by leading equipment manufacturers

Active ingredients: Sodium percarbonate peroxohydrate Available as: 500 g bag Our recommendation: OroCup, MD 555 cleaner, MD 550

- Dissolves encrustations of lime scale, prophy powder or prophylaxis pearl products
- Highly intensive cleaning effect against deposits and sludge
- Certified material compatibility recommended by leading equipment manufacturers
- Ideal in combination with Orotol<sup>®</sup> plus

Ingredients: Organic and inorganic acids, foam-free surfactants Available as: 2.5 Our recommendation: OroCup, Orotol® plus, Orotol® ultra, MD 550

agents

- Simultaneous cleaning, deodorisation and care of spittoons
- No abrasive components the spittoon suffers no scratching
- Anti-microbial effective against bacteria and yeast fungus











The origina











# Impressions and dentistry work

## MD 530 Cement remover and denture cleaner

Ready-to-use solution for removal of cement, calculus, smoking stains, oxide layers etc. from dentures, bridges, crowns, spatulas, glass etc.



## MD 535 Plaster remover

Ready-to-use solution for removal of plaster and alginates from impression trays and other aids



- Intensive cleaning effect and residue-free removal of acid-soluble plaque
- For direct use, outstanding material compatibility
- Particularly fast and effective, thanks to a special formulation
- Contains corrosion inhibitors
- Also suitable for use in an ultrasonic bath

Ingredients: Inorganic acids, surfactants Available as: 2.5 | Our recommendation: Hygobox, Hygosonic

- Neutral special cleaner with particularly fast action
- Intensive cleaning power for residue-free removal of plaster and alginate
- With activity indicator
- Also suitable for use in an ultrasonic bath

Ingredients: Alkaline cleaning components and complexing agents Available as: 2.5 | Our recommendation: Hygobox, Hygosonic

#### MD 520 Impression disinfection

Alcohol-free ready-to-use solution for effective and gentle disinfection and cleaning of impressions, prosthetic work, etc.



- Bactericidal (incl. TB) yeasticidal, virucidal (certified to EN 14476)

- For use with alginates, silicones, polyether rubber, hydrocolloids, polysulphides etc.
- Increased accuracy in the plaster master model, thanks to effective cleaning effect
- No adverse effects on dimensional stability or plaster compatibility

Active ingredients: Quaternary ammonium compounds, glutaraldehyde Available as: 2.5 |, 10 | Our recommendation: Disinfection box, Hygobox, Hygojet

# Unit

### Hygojet

Disinfection and cleaning sheath for impressions and prosthetic work



- The closed system prevents cross-infections
- Economical to use with just 10 ml MD 520 required per impression
- Extremely effective for 6 impressions in 12 minutes
- Tested disinfection process proven by experts

Dimensions (WxHxD): 56 x 34 x 32 cm Our recommendation: MD 520

# Spectrums of activity and applications

	Disinfection	Cleaning	Aldehyde-free	Biodegradable	VAH listed	IHO list of disinfectants	Tested in accordance with EU standards	4-field test in accordance with VAH/EN 16615	Bactericidal	Tuberculocidal	Yeasticidal	Fungicidal	Limited virucidal	Limited virucidal plus	Virucidal	Application	Action time in minutes as per VAH	Reaction time in minutes in an ultrasonic bath	No. of days product works for
Instruments																			
ID 212 Instrument disinfection									<b>5</b>	<b>5</b>	<b>5</b>	i				2%	5	5 <sup>10</sup>	7
ID 212 forte Instrument disinfection									<b>5</b>	<b>5</b>	<b>=</b> <sup>5</sup>	;				2%	5	5 <sup>8</sup>	7
ID 213 Instrument disinfection									<b>5</b>	<b>5</b>	<b>5</b>	∎ <sup>5</sup>			■ <sup>16</sup>	2%	5/60 <sup>7</sup>	1°	14
ID 215 Enzymatic instrument cleaner																1%/2%	30/15	15/1	
ID 220 Bur disinfection									<b>5</b>	<b>5</b>	<b>5</b>	∎ <sup>5</sup>			■ <sup>15,16</sup>	Undiluted	17	17	7
Surfaces																			
FD 300 Surface disinfection								∎ <sup>5</sup>	<b>5</b>	<b>5</b>	<b>5</b>	i			■ <sup>15,16</sup>	1%/2%	1510/5		
FD 312 Surface disinfection <sup>1</sup>								∎ <sup>5</sup>	<b>5</b>		<b>5</b>	i				1%	15		
FD 370 cleaner Cleaning agent for medical practices																1%/2%			
FD 322 Quick-acting disinfection								5,6	<b>■</b> <sup>6</sup>	<b>6</b>		<sup>6</sup> ∎				Undiluted	17		
FD 333 Quick-acting disinfection								<b>■</b> <sup>6</sup>	<b>■</b> <sup>6</sup>	<b>6</b>	<b>•</b>	<sup>6</sup> ∎				Undiluted	17		
FD 333 forte Quick-acting disinfection								■ <sup>5,6</sup>	<b>■</b> <sup>5,6</sup>	<sup>5</sup> ∎ <sup>5,6</sup>	5,0 ■ <sup>5,0</sup>	<sup>6</sup> <b>=</b> <sup>5,6</sup>			15,16	Undiluted	17		
FD 366 sensitive Disinfection of sensitive surfaces				-				6	<b>■</b> <sup>6</sup>	<b>■</b> <sup>6</sup>	∎ <sup>6</sup>		-			Undiluted	րո		
FD 300 top wipes Surface disinfection					■ <sup>17</sup>			<b>5</b>	<b>5</b>	<b>5</b>	<b>5</b>	<b>∎</b> <sup>5</sup>			<b>■</b> <sup>16</sup>	Ready-to- use	2 <sup>10</sup>		<u> </u>
FD 312 wet wipes Surface disinfection <sup>1</sup>								<b>5</b>	<b>5</b>		<b>5</b>	i				Ready-to- use	15		<u> </u>
FD 322 premium wipes quick-acting disinfection / FD 322 top wipes quick-acting disinfection					-		-	5,6	∎ <sup>5,6</sup>	5,6	• ■ <sup>5,0</sup>	<sup>5</sup> ∎ <sup>5,6</sup>				Ready-to- use	17		
FD 333 wipes Quick-acting disinfection FD 333 top wipes Quick-acting disinfection	•				■ <sup>17</sup>	■ <sup>17</sup>	-		<b>■</b> <sup>6</sup>	<b>•</b>	<b>■</b> <sup>6</sup>	<b>6</b>	•	•		Ready-to- use	17		
FD 333 forte premium wipes quick-acting disinfection, FD 333 forte wipes quick-acting disinfection	-				■ <sup>17</sup>		-	5,6	<b>■</b> <sup>5,6</sup>	<b>■</b> <sup>5,6</sup>	∎ <sup>5,0</sup>	<sup>5,6</sup>		-	■ <sup>15,16</sup>	Ready-to- use	17		
FD 360 Cleaning and care of vinyl upholstery																Ready-to- use			
FD 350 Disinfection wipes									<b>■</b> <sup>6</sup>	<b>6</b>		• <b>■</b> <sup>6</sup>				Ready-to- use	17		
FD 366 sensitive wipes Disinfection of sensitive surfaces	•						•	5,6	<b>•</b> 6	<b>■</b> <sup>6</sup>	•	i				Ready-to- use	1 <sup>11</sup>		
FD 366 sensitive top wipes Disinfection of sensitive surfaces					•	■ <sup>17</sup>	•	5,6	∎ <sup>5,6</sup>	∎ <sup>5,6</sup>	∎ <sup>5,0</sup>	ò				Ready-to- use	1 <sup>11</sup>		
Skin and hands																			
HD 410 Hand disinfection <sup>1, 2, 3</sup>									<b>•</b> <sup>6</sup>	<b>6</b>	<b>6</b>					Undiluted	0,5/1,512,13	1	
HD 425 Cleansing lotion <sup>4</sup>																Undiluted			
HD 435 Cleansing lotion <sup>4</sup>																Undiluted			
HD 440 Care lotion																Undiluted			
Special areas																			
Orotol <sup>®</sup> plus disinfection of suction systems								∎ <sup>5</sup>	<b>5</b>		<b>∎</b> <sup>5</sup>	5				2%	5/60		
Orotol® ultra disinfection of suction systems									<b>5</b>	<b>5</b>	<b>5</b>	i ∎ <sup>5</sup>			<b>1</b> 5,16	1%	5		
MD 555 cleaner Special detergent for suction systems																5%	30-120		
MD 550 Spittoon bowl cleaner	+								<b>6</b>		<b>■</b> <sup>6</sup>	,				Undiluted	approx. 5		<u> </u>
MD 520 Impression disinfection	+								<b>5</b>	6	5	;			15,16	Undiluted	5/1014		7
-					_	_									_	onanoioa	3/10		
MD 530 Cement remover / denture cleaner	+-			-	-	-	-		-	-	-		-		_	Undiluted	Max. 10	Max. 10	

<sup>1</sup>Use biocides safely. Always read the label and product information before use. <sup>2</sup>Lipid replenishing <sup>3</sup>RKHisted <sup>4</sup>pH value typical for skin <sup>5</sup>High burden<sup>6</sup>Low burden <sup>7</sup>Incl. TB bacteria <sup>8</sup>For TB bacteria 4%/30 min. <sup>9</sup>For TB bacteria 2%/30 min.<sup>10</sup>For TB bacteria 60 min. <sup>11</sup>For TB bacteria 2 min. <sup>12</sup>For TB bacteria 30 sec. <sup>13</sup>Hygienic hand disinfection/surgical hand disinfection <sup>14</sup>In the Hygojet <sup>15</sup>Certified to DVV/RKI <sup>16</sup>Certified to DVV/RKI <sup>16</sup>Certified to EN 14476 <sup>17</sup>Active ingredient solution

DÜRR DENTAL SE Höpfigheimer Str. 17 74321 Bietigheim-Bissingen Germany www.duerrdental.com info@duerrdental.com

