



www.rimsa.it

Brightening ideas

Designers of the world's first scialytic lamp using LED technology, RIMSA now applies the experience it has gained over the years to its new **U29** Model.

U29 combines the indirect technology consolidated by the success obtained by the Pentaled Series with an elegant and compact design. U29 is a high-performance lamp, strongly engineered but accessible, so that any surgeon can benefit from it.

The U29 lamp has been developed to provide a comfortable light without any glare for surgeon and assistants, with an extremely easy-to-handle dome and a reasonable price. This product therefore goes to join what is an already broad range of RIMSA surgical lamps.





Controls The controls on the

- membrane keyboard activate the following functions:
- colour temperature selection (5000°K or 4500°K)
- light intensity adjustment (up to 160 Klx)
- light field adjustment
 switching on/off
- EndoLed, a special type of light that facilitates endoscopic operations. The surgeon can select the exact brightness and colour temperature required.

Dome

- The dome, which is very compact and has a large light emission surface, is balanced so as to make its operation light and stable.
 It features an easy-to-clean shape and is without screws and surface holes.
- The dome features an ergonomic round grip to hold and move the lamp.
- The screen roughness collimates the light beam and improves the scialytic effect.
 The polycarbonate screen is treated to be resistant to scratches and knocks.

indirect light



Indirect light, a crucial feature of the U29 Model, provides the surgeon with cold and deep light, without any sort of glare. This is due to especially calculated parabolic dishes that reflect the entire spectrum of light emitted by each LED lamp in a punctiform manner, without dispersion. This means:

- fewer LEDs needed
- lower irradiated heat in the operating field

Thanks to the physical principle of indirect light, the surgeon and his/her assistants are never disturbed by the light emitted and never become dazzled. The light is always on focus at different distances (from 80 cm to 150 cm). The size of the light field can be adjusted electronically according to the type of light required by the surgical speciality.



U29 mobile	U29+ 29+DY2	U29 ceiling			
	RIASA RIASA RIASA	RIMSA o			

Performances:

PERFORMANCES U29		ELECTRICAL DATA	
Light intensity at 1m distance (E _c)	160 Klx	Primary alternating voltage (V a.c.)	100 ÷ 240
Color temperature: double selection	4500 - 5000 K	Secondary continue voltage (V d.c)	24
Color rendering index (CRI)	96 Ra	Frequency (Hz)	50/60
Color rendering index (R9)	≥ 90	Absorbed power	48 W - 56 VA
No. of Leds	29		
Light field diameter adjustment	Electronic	GENERAL DATA	
Light field diameter adjustable from – to –	140 – 230 mm	Colour	RAL 9003
d10 light field diameter where illuminance reaches 10% of Ec	200 mm	Directive	93/42/EEC*
d50 light field diameter where illuminance reaches 50% of Ec	120 mm	Standards Classification of product (Medical Device)	IEC60601-2-41 Class I
Depth of illumination IEC 60601-2-41 (L1+L2) at 60%	490 mm	DIMENSIONS	
Depth of illumination IEC 60601-2-41 (L1+L2) at 20%	1350 mm	Out reflector diameter (cm)	52
Control of the illuminance	20 - 100 %	Useful lighting surface (cm2)	737
Total radiated energy Ee where the illuminance reaches max level (W/m²)	580	OPTIONAL	
Ratio between radiated energy Ee and illuminance Ec (mW/m²lx)	3,68	Battery group with automatic charger	Optional
Radiated UV energy with wavelength less than 400nm (W/m²)	0,001	HD camera	Optional
		IR remote control	Optional
Average LED life	> 60.000 hours	Wall control panel	Optional
		Monitor arm	Optional

All lighting values are subjected to a tolerance of \pm 5% due to manufacturing and metrological reasons. * main directive and further amendments



RIMSA P. LONGONI SRL

Via Monterosa, 18/22 20831 Seregno (MB) - Italy Tel. + 39 0362 325709 Fax + 39 0362 328559 E-mail: info@rimsa.it Rimsa retain a right to improve the products in the catalogue without notice. Reproduction in part or in whole of any product is forbidden.



www.rimsa.it

