

I LIKE IT!

Surfaces – sparkling clean with maximum safety



The thorough but gentle cleaning and disinfection of surfaces eliminates risks while helping to ensure that the surgery's equipment retains its value. Dürr Dental offers a whole host of concentrates and ready-to-use solutions to assist you with this, along with an extensive range of wipes.

Alcohol-free disinfection and cleaning



FD 300 | Surface disinfection

Fully virucidal and alcohol-free concentrate for cleaning and disinfecting surfaces of non-invasive medical devices and medical fixtures

- Especially fast reaction time
- Economical application: Up to 250 l of ready-to-use solution at 1 %
- Good material compatibility

Spectrum of activity¹

Bactericidal, yeasticidal, tuberculocidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), virucidal (against enveloped viruses and non-enveloped viruses)

Active agents

Alkylamines, quaternary ammonium compounds

Units

2.5 litre

Ideal addition

Dry wipes (see p. 20)



FD 312 | Surface disinfection*

Alcohol-free concentrate for cleaning and disinfecting surfaces of non-invasive medical devices, fixtures and floors

- Economical application: Up to 250 l of ready-to-use solution at 1 %
- Good material compatibility
- Pleasant scent

Spectrum of activity²

Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses)

Active agents

Alkyl(C12-18)dimethylbenzylammonium chloride (quaternary ammonium compounds)

Units

2.5 litre

Ideal addition

Dry wipes (see p.20)



FD 312 wet wipes | Surface disinfection*

Alcohol-free disinfection and cleaning wipes for floors, surgery fixtures, etc.

- High-yield surface coverage
- Wipes with tear-resistant honeycomb structure
- Pleasant scent

Spectrum of activity²

Bactericidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses)

Active agents

Alkyl(C12-18)dimethylbenzylammonium chloride (quaternary ammonium compounds)

Units

Package with 15 wipes

Wipe size

22 x 42 cm



¹ Tested in accordance with EN 13727, EN 13624, EN 14348, EN 14476 and EN 16615.

² Tested in accordance with EN 13727, EN 13624, EN 14476 and EN 16615.

* Use disinfectants safely. Always read the label and product information before use.

Alcohol-free ready-to-use solutions for disinfection



Perfume-free option also available

FD 322 | Quick-acting disinfection

Fast-acting, ready-to-use solution for the disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Very short reaction time
- Dries without residue
- Good material compatibility

Spectrum of activity¹

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40

Active agents

Alcohols

Units

750 ml for Hygowipe Plus, 1 litre, 2.5 litres, 10 litres

Ideal addition

FD 322 perfume-free: 2.5 litres, 10 litres

Dry wipes (see p. 20), 600 ml spray bottle, spray head



FD 333 | Quick-acting disinfection

Fast-acting, ready-to-use solution for disinfecting alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Good material compatibility
- Fresh lemon scent
- Quick drying

Spectrum of activity¹

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses)

Active agents

Alcohols, quaternary ammonium compounds

Units

750 ml for Hygowipe Plus, 1 litre, 2.5 litres, 10 litres

Ideal addition

Dry wipes (see p. 20), 600 ml spray bottle, spray head



FD 333 forte | Quick-acting disinfection

Fully virucidal ready-to-use solution for the disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Good material compatibility
- Fresh lemon scent
- Quick-acting disinfection effect

Spectrum of activity¹

Bactericidal, tuberculocidal, mycobactericidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses) virucidal (enveloped and non-enveloped viruses)

Active agents

Alcohols, quaternary ammonium compounds

Units

1 litre, 2.5 litres, 10 litres

Ideal addition

Dry wipes (see p. 20), 600 ml spray bottle, spray head



Reduced alcohol content

FD 366 sensitive | Disinfection of sensitive surfaces

Alcohol-reduced ready-to-use solution for the quick disinfection of sensitive surfaces of invasive and non-invasive medical devices and fixtures

- Quick-acting disinfection effect
- Certified material compatibility – recommended by leading equipment manufacturers

Spectrum of activity¹

Bactericidal, tuberculocidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses

Active agents

Alcohols, quaternary ammonium compounds

Units

750 ml for Hygowipe Plus, 1 litre, 2.5 litres, 10 litres

Ideal addition

Dry wipes (see p. 20), 600 ml spray bottle, spray head, FD 360 Cleaning and care of vinyl upholstery



Alcohol-based disinfection wipes

For alcohol-resistant surfaces, available in 3 sizes



FD 322 wipes | Quick-acting disinfection

Disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Very short reaction time
- Dries without residue
- Good material compatibility

Spectrum of activity ¹

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40

Active agents Units

Alcohols
Dispenser box of 100 wipes,
refill pack of 100 wipes

Wipe size

20 x 30 cm



FD 322 top wipes | Quick-acting disinfection

Product details, see above

Units

Package with 50 wipes

Wipe size

17.5 x 25 cm



FD 322 premium wipes | Quick-acting disinfection

Product details, see above

Units

Dispenser box of 80 wipes,
refill pack of 80 wipes

Wipe size

14 x 19 cm

Material compatibility of disinfectants

A disinfectant with a high alcohol content can damage sensitive materials and impair their durability. Therefore, use a suitable disinfectant for each type of material.

More info in the video at duerrdental.com/FD366



Alcohol-based disinfection wipes

Perfume-free, available in 2 sizes



Perfume-free

FD 333 wipes | Quick-acting disinfection

Disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Quick-acting disinfection effect
- Quick drying times
- Perfume-free

Spectrum of activity¹

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses)

Active agents
Units

Alcohols, quaternary ammonium compounds
Dispenser box of 110 wipes, refill pack of 110 wipes

Wipe size

14 x 22 cm



Perfume-free

FD 333 top wipes | Quick-acting disinfection

Product details, see above

Units

Package of 50 wipes

Wipe size

17.5 x 25 cm

Fully virucidal, available in 2 sizes



FD 333 forte wipes | Quick-acting disinfection

Fully virucidal, ready-to-use disinfection wipes for fast disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and inventory

- Quick-acting disinfection effect
- Fresh lemon scent

Spectrum of activity¹

Bactericidal, tuberculocidal, mycobactericidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), virucidal (enveloped and non-enveloped viruses)

Active agents
Units

Alcohols, quaternary ammonium compounds
Dispenser box of 100 wipes, refill pack of 100 wipes

Wipe size

20 x 30 cm



FD 333 forte premium wipes | Quick-acting disinfection

Product details, see above

Units

Dispenser box of 80 wipes, refill pack of 80 wipes

Wipe size

14 x 19 cm



For small areas, available in 2 variants



FD 350 | Quick-acting disinfection



Disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Quick-acting disinfection effect
- Quick, streak-free drying
- 3 different fragrances:
- Classic, Flower, Lemon

Spectrum of activity¹

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40

Active agents

Alcohols

Units

Dispenser box of 110 wipes, refill pack of 110 wipes

Wipe size

14 x 22 cm



Plastic-free natural fibres

FD 350 green | Quick-acting disinfection



Ready-to-use disinfection wipes made from plastic-free natural fibres for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Quick-acting disinfection effect
- Quick, streak-free drying
- Good material compatibility and particularly large coverage per wipe

Spectrum of activity¹

Bactericidal, tuberculocidal, fungicidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), limited virucidal PLUS (incl. rotaviruses, adenoviruses, noroviruses), polyomaviruses SV40

Active agents

Combination of alcohols

Units

Dispenser box of 50 wipes, refill pack of 50 wipes

Wipe size

14 x 19 cm

For alcohol-resistant surfaces, available in 2 sizes



Reduced alcohol content

FD 366 sensitive wipes | Quick-acting disinfection

Ready-to-use disinfection wipes for the quick disinfection of alcohol-resistant surfaces of invasive and non-invasive medical devices and fixtures

- Quick-acting disinfection effect
- Recommended by equipment manufacturers
- Perfume-free

Spectrum of activity¹

Bactericidal, tuberculocidal, yeasticidal, limited virucidal (enveloped viruses incl. HBV, HCV, HIV and coronaviruses), noroviruses

Active agents

Alcohols, quaternary ammonium compounds

Units

Dispenser box of 100 wipes, refill pack of 100 wipes

Wipe size

20 x 30 cm



Reduced alcohol content

FD 366 sensitive top wipes | Quick-acting disinfection

Product details, see above

Units

Package with 50 wipes

Wipe size

17.5 x 25 cm

Dry wipes

Ideal for large surfaces, available in 2 sizes



Suitable for universal use

FD multi wipes

Universal dry wipes, can be soaked in disinfectant and cleaning agents

- All-purpose wipes can be used together with surface disinfectants
- Lint-free and tear-proof
- Excellent wetting; does not scratch surfaces

Units

Dispenser box of 180 wipes, refill pack of 180 wipes

Wipe size

20 x 30 cm

Ideal addition

Can be used for surface disinfection with all Dürre Dental ready-to-use solutions and concentrates. Tested and approved using the 4-field test (EN 16615).



Suitable for universal use

FD multi wipes compact

Product details, see above

Units

Dispenser box of 60 wipes, refill pack of 60 wipes

Wipe size

29 x 30 cm

Ideal addition

Can be used for surface disinfection with all Dürre Dental ready-to-use solutions and concentrates. Tested and approved using the 4-field test (EN 16615).

Eco-friendly wipes, available in 2 sizes



Plastic-free natural fibres

FD multi wipes green

Universal dry wipes for surface disinfection made from plastic-free natural fibres using eco-friendly production methods

- Very economical with efficient coverage
- Tear-resistant, lint-free and finely structured non-woven wipes
- Excellent wetting; does not scratch surfaces

Units

Dispenser box of 180 wipes, refill pack of 180 wipes

Wipe size

20 x 30 cm

Ideal addition

FD 312 (4%), FD 322, FD 333, FD 333 forte. Tested and approved using the 4-field test (EN 16615).



Plastic-free natural fibres

FD multi wipes compact green

Product details, see above

Units

Dispenser box of 60 wipes, refill pack of 60 wipes

Wipe size

29 x 30 cm

Ideal addition

FD 312 (4%), FD 322, FD 333, FD 333 forte. Tested and approved using the 4-field test (EN 16615).



Cleaning and care of vinyl upholstery

Recommended by equipment manufacturers



FD 360 Cleaning and care of vinyl upholstery

Cleaning fluid for deep cleaning and care of vinyl upholstery

- Cleaning agent, gentle on materials, combined with a care product for long-term value-retention
- Intensive cleaning effect – removes stains of all kinds, rims and discolouration
- Certified material compatibility – recommended by leading equipment manufacturers
- Ideal for treatment units in combination with FD 366 sensitive

Ingredients

Special surfactants, silicon compounds and care formula based on avocado oil

Units

Set with 250 ml spray bottle, 1 FD polishing wipe and 5 special sponges

Ideal addition

FD 366 sensitive, FD 366 sensitive wipes, FD 366 sensitive top wipes, FD polishing wipe



FD polishing wipe

Microfibre cloth with special fibres for gentle cleaning and polishing of sensitive surfaces

- Particularly thorough thanks to the fine fibre material
- Fast, streak-free cleaning
- Antistatic properties remove small particles of dirt and dust
- Suitable for washing at temperatures of up to 60 °C

Units

Available as individual wipes

Wipe size

32 x 32 cm

Our products for eco-friendly surface disinfection

The pre-soaked FD 350 green Disinfection wipes as well as the dry non-woven wipes FD multi wipes green and FD multi wipes compact green are made from plastic-free natural fibres using eco-friendly production methods.

green



Scan the QR code for more information.

